


MICROSCOPES/YEAR SOLD IN EUROPE



Biological Microscope Series

The perfect entry into the world of microscopy; robust, easy to handle, and affordable: The Motic SFC100 Series offers the ideal instrument for a first step into the "world of small things". No plastics, but efficient optics 100% made of coated glass, this is the way to crisp and clear

images for a long-term usage in classroom or for leisure activities. The SFC100 is a classic microscope with a classic performance.



ST30C 2L00

ST30C 6LED

ST36C 2L00

ST36C 6LED

Stereo Microscope Series

This microscope series contains 4 models. They are developed for a first contact with tiny objects, not visible to the unarmed human eye, but ready to be shown to a younger public. A rear-facing head model, a new forward-facing head model: depending on your choice both

constructions are available. The cordless LED illumination will probably rule the standard tungsten illumination, as its ecological advantages are obvious: less heat development, less power consumption.

SFC100 Series Specifications	SFC100 FL	SFC100 FLED
Observation Tube		
45° Inclined Monocular head for comfortable viewing while seated, rotates 360°	•	•
Eyepieces		
Widefield eyepiece 10X/18mm with built-in pointer	•	
Nosepiece		
Revolving triple nosepiece with positive click-stops for precise magnification changes	•	•
Objectives		
DIN Achromatic objectives parfocal and parcentric 4X/0.1, 10X/0.25 and 40X/0.65 with spring	•	•
Focus mechanism		
Separate coarse and fine focusing adjustment with upper-limit stop to prevent damage to slides and objectives	•	•
Stage		
Large 120mmx110mm stage with locked-on, spring stage clips	•	•
Condenser		
Special built-in 0.65NA condenser with iris diaphragm and filter holder		•
Illumination		
Built-in tungsten illumination	•	
Rechargeable LED illumination, 3.5V 70mW with intensity control, battery autonomy of 50 hours	T BOOK TO SEE	•
MSRP Manufacturer Suggested Retail Price	134,00€	157,00€

ST30 Series Specifications	ST30C 2L00	ST30C 6LED	ST36C 2L00	ST36C 6LED
Observation Tube	Assistante	E CONTRACTOR		VAROUS SEA
45° Inclined rear-facing binocular head with diopter adjustment on left eyetube	•	•		
45° Inclined 360° rotatable binocular head with diopter adjustment on left eyetube	3000		•	•
Eyepieces		TO SHE		
Widefield eyepieces 10X/20mm	•	•	•	•
Objectives		NUMBER OF		
2X/4X magnification changer. Objectives are parfocal, parcentric and achromatic	•	•	•	• /
Focus mechanism				
Heavy-duty rack & pinion focusing block with slip clutch	•		•	•
Stage Plates		AND TO B		
Supplied with 80mm stage plates (frosted glass and black/white)	•	•	•	•
Illumination		SAME TO BE		
Pole-stand with incident and transmitted 12V/10W tungsten illumination	•		•	AT STATE OF
Pole-stand with incident and transmitted rechargeable LED illumination 3.1V 60mW per piece, battery autonomy of 50 hours	1939/2			1
MSRP Manufacturer Suggested Retail Price	151,00€	189,00€	165,00€	205,00€

The prices shown are recommended retail prices and may vary depending on the country. Taxes not included.

ST30 Series Accessories

SFC100 Series Accessories

Eyepieces

Huygens eyepiece H5X / 14.5 mm Widefield eyepiece WF10X / 18 mm Widefield eyepiece WF15X / 12 mm Widefield eyepiece WF20X / 11 mm

Micrometer eyepieces

Micrometer eyepiece WF10X / 18mm, 100 divisions in 10mm

Achromatic Objectives

Achromatic objective 4X / 0.10
Achromatic objective 10X / 0.25
Achromatic objective 20X / 0.40
Achromatic objective 40X / 0.65 / S
Achromatic objective 60X / 0.85 / S
Achromatic objective 100X / 1.25 / S-oil

Polarizing equipment

Polarizing equipment B.201

Spare parts

Mechanical stage MSH-009
Blue filter
Dust cover (size B)
42 mm plano-concave mirror and fork
Immersion oil 5ml
Plastic carrying case
Aluminium box
Cleaning kit
Lens tissue paper (Pack 500)
Motic slides 76x25mm (50 units)
Motic coverslips 22x22mm (100 units)

Spare bulbs

220V/20W Tungsten bulb LED 20mA, 3.5V, 70mW

Eyepieces

Widefield eyepiece WF5X / 22 mm Widefield eyepiece WF10X / 20 mm Widefield eyepiece WF15X / 13 mm Widefield eyepiece WF20X / 10 mm Widefield eyepiece WF30X / 8 mm

Micrometer eyepieces Micrometer eyepiece WF10X/20mm,

140 divisions in 14mm and crosshair Micrometer eyepiece WF10X/20mm, 70 divisions in 14mm Micrometer eyepiece WF10X/20mm, 360° protractor with 1° divisions and crosshair Micrometer eyepiece WF10X/20mm, 100 divisions in 10mm

Darkfield and Polarizing equipment

Darkfield attachment with iris diaphragm Jewellery clamp Polarizing equipment B.301 (1-piece rotary)

Spare parts

White stage plate, Ø 80 mm
Black stage plate, Ø 80 mm
Frosted glass stage plate, Ø 80 mm
Dust cover (size B)
Eyecup (unit)
Plastic carrying case
Aluminium box
Cleaning kit
Lens tissue paper (Pack 500)

Spare bulbs

12V/10W Tungsten bulb LED incident light assembly LED transmitted light assembly MAKE YOUR LESSONS
MORE FUN

DETICAL WARRANTY

MORETHAN

SFC100 ERIES ST30 SERIES MICROSCOPESNEAR!

SOLD IN EUROPE K





Canada | China | Germany | Spain | USA



www.moticeurope.com

WETZLAR, GERMANY

SALES OFFICE

Christian-Kremp-Straße 11, 35578 Wetzlar, Germany T. +49 6441 21001 0 | Fax: 49-6441-210 0122

BARCELONA, SPAIN

SALES OFFICE & EUROPEAN DISTRIBUTION CENTRE

C. Les Corts 12, Pol. Ind. Les Corts 08349 Cabrera de Mar, Barcelona, Spain T. +34 93 756 62 86 | Fax: 34 93 756 62 87







